


What do trend arrows mean?

Expected blood sugar variation within the next 10 minutes according to the trend arrows.

	Dexcom G6 G7	Medtronic Guardian Connect Guardian 3 Guardian 4	Abbott FreeStyle Libre 2
		Rapid increase higher than 1.7 mmol/L	
	Rapid increase higher than 1.7 mmol/L	Increase 1.1 - 1.7 mmol/L	
	Increase 1.1 - 1.7 mmol/L	Slow increase 0.5 - 1.0 mmol/L	Rapid increase higher than 1.0 mmol/L
	Slow increase 0.6 - 1.1 mmol/L		Increase 0.6 - 1.0 mmol/L
 or no arrow	Stable increase or decrease lower than 0.6 mmol/L	Stable increase or decrease lower than 0.6 mmol/L	Stable increase or decrease lower than 0.6 mmol/L
	Slow decrease 0.6 - 1.1 mmol/L		Decrease 0.6 - 1.0 mmol/L
	Decrease 1.1 - 1.7 mmol/L	Slow decrease 0.5 - 1.0 mmol/L	Rapid Decrease higher than 1.0 mmol/L
	Rapid Decrease higher than 1.7 mmol/L	Decrease 1.1 - 1.7 mmol/L	
		Rapid Decrease higher than 1.7 mmol/L	